TRACOE medical GmbH  
Reichelsheimer Straße 1 / 3, 55268 Nieder-Olm / Germany  
Tel.: +49 6136 9169-0, Fax: +49 6136 9169-200  
info@tracoe.com, www.tracoe.com

**TRACOE**<sup>®</sup>  
*technic*

**CE**  
0 1 2 3  
2018-07  
10456  
90

**TRACOE<sup>®</sup> Connection Tubes for TRACOE<sup>®</sup> cuff pressure monitor (REF 720)/ TRACOE<sup>®</sup> cuff pressure monitor sensitive (REF 721) and TRACOE<sup>®</sup> smart Cuff Manager (REF 730)**

**TRACOE<sup>®</sup> Verbindungsschläuche für TRACOE<sup>®</sup> cuff pressure monitor (REF 720)/ TRACOE<sup>®</sup> cuff pressure monitor sensitive (REF 721) und TRACOE<sup>®</sup> smart Cuffmanager (REF 730)**

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**EN**

**Instructions for Use for TRACOE<sup>®</sup> Connection Tubes for TRACOE<sup>®</sup> cuff pressure monitor (REF 720)/ TRACOE<sup>®</sup> cuff pressure monitor sensitive (REF 721) and TRACOE<sup>®</sup> smart Cuff Manager (REF 730)**

The following instructions are for the information of the doctor and nursing staff. They relate exclusively to the TRACOE connection tubes for TRACOE cuff pressure monitor/ TRACOE cuff pressure monitor sensitive as well as for the TRACOE smart Cuff Manager (REF 702 and REF 707) and must be read carefully before use.

#### 1. Pack Contents

The package of the REF 702 contain 10 units and that of the REF 707 contains 5 units of individually packed sterile connection tubes for the TRACOE cuff pressure monitor/ TRACOE cuff pressure monitor sensitive or TRACOE smart Cuff Manager.

All tubes are fitted with a female Luer lock connector with check valve at one end and a male Luer lock connector at the other end. Lengths of the connection tubes:  
REF 702: 1.0 m  
REF 707: 0.3 m

#### 2. Intended Use

The connection tubes are designed for the pneumatic connection between the TRACOE cpm (REF 720)/ TRACOE cpm sensitive (REF 721) and TRACOE smart Cuff Manager (REF 730) respectively and the inflation tube of a tracheostomy tube or an endotracheal tube.

**Single-patient Product:** REF 702 and REF 707 are medical products intended for only one patient. Multiple use on the same patient is allowed.

**Useful life:** The maximum useful life is 72 hours. Products that show only minor damage must be replaced.

#### 3. Sterilisation

The connection tube is sterilised with ethylene oxide. Sterility is guaranteed only if the pack is intact and sealed all round. The connection tube must not be resterilised after use.

#### 4. Instructions

To create a pneumatic connection between the TRACOE cuff pressure monitor/ TRACOE cuff pressure monitor sensitive or the

TRACOE smart Cuff Manager and the cuff of a tracheostomy tube or an endotracheal tube, the connection tube with its female Luer lock connector with check valve is pushed firmly onto the corresponding male component on the TRACOE cuff pressure monitor/ TRACOE cuff pressure monitor sensitive or the TRACOE smart Cuff Manager. The male Luer lock connector located on the other end of the tube is pushed or screwed firmly onto the valve of the inflation line supplying the cuff on the tracheostomy tube or endotracheal tube. The seal integrity of the connection is then tested using the TRACOE cuff pressure monitor/ TRACOE cuff pressure monitor sensitive.

#### 5. Returns

Returns of used products are accepted only if accompanied by a completed decontamination certificate. You can obtain the required form either directly from TRACOE medical or from the website: www.tracoe.com.

#### 6. General Terms and Conditions

The sale, delivery and return of all TRACOE products are carried out exclusively according to the General Terms and Conditions which can be obtained directly from TRACOE medical GmbH or from the website: www.tracoe.com.

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