cuff pressure monitor sensitive

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English

INTENDED USE

The Cuff Pressure Gauge is used for inflation, deflation and pressure control of tracheal tubes with high volume low pressure cuffs. The vacuum valve is used for deflation of low pressure cuffs.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by trained medical personnel.
- Before use, the products must be inspected visually for damages (cracks, breakage, etc.) and a functional check must be performed (see section "Functional check").
- The cuff pressure gauge is not suitable for MRI!
- In case that the cuff pressure gauge is used for monitoring and remains on the inflation line, the pressure must be checked regularly every 15 minutes by the user. Pressure increase can be adjusted by pressing the red release valve and pressure decrease can be adjusted by squeezing the inflation bulb. This cuff pressure gauge is a manual system and can not hold the cuff pressure constantly over a longer period. For such case the manufacturer offers an automatic Cuff Controller which monitors the cuff pressure electronically and maintains the pre-set pressure.
- When disconnecting the cuff pressure gauge from the inflation line, it has to be disconnected together with the connecting tube as otherwise the system is open.
- A mechanical stop is integrated in the housing to prevent that the release valve gets stuck. The release valve can get stuck by excessive strong activation. In this case, the device must not be used further and must be returned to the manufacturer for repair.
- The connecting tube is for single use and must not be reprocessed. Reprocessing influences the function of the connecting tube. Any reuse involves the potential risk for cross contamination.

FUNCTIONAL CHECK

Test the Cuff Pressure Gauge as follows:

1a) Seal Luer nozzle with the finger (picture 1).

1b) Inflate with the inflation bulb to 40 cmH₂O; the value must be constant for 2 - 3 sec.; in case of pressure decrease, the device must be returned to the manufacturer for repair.

If a connecting tube is used, the system consisting of Cuff Pressure Gauge and connecting tube has to be checked (additionally to 1):

2a) Attach the connecting tube to the Luer nozzle of the Cuff Pressure Gauge (picture 2).

2b) Seal the end of the connecting tube with the finger (picture 3).

2c) Inflate with the hand inflation bulb to 40 cmH₂O; the value must be constant for 2 - 3 sec.; in case of a pressure decrease, the connecting tube is leaking and must be replaced. Check again with the new connecting tube (repeat steps 2a - 2c).

During connection and disconnection of the connecting tube, do not bend the luer fitting of the device, to prevent damage. Furthermore, make sure that the connecting tube is held on the connector (not at the tubing).

If the vacuum valve is used, do the following test prior to use:

- 3a) Inflate the low pressure cuff of a tracheal tube with the Cuff Pressure Gauge to 40 cmH₂O (picture 4).
- 3b) Attach the inflation line to the Luer nozzle of the red vacuum valve (picture 5, 6).
- 3c) Squeeze the inflation bulb several times; the low pressure cuff has to deflate.
- 3d) Squeeze the inflation bulb until it does not return to its original position; the low pressure cuff must be completely empty now (picture 6).
- 3e) Disconnect the inflation line from the vacuum valve. Connect the inflation line to the Luer nozzle of the Cuff Pressure Gauge. The needle of the Cuff Pressure Gauge has to show a negative pressure by moving counter-clockwise (picture 7).

INSTRUCTION FOR USE

Before use check the low pressure cuff for leaks. Prior to intubation or extubation withdraw all air from the cuff via a syringe or the vacuum valve. Connect the Cuff Pressure Gauge (eventually by using the connecting tubing) to the inflation line of the tracheal tube. Inflate the cuff to a pressure of 60 - 90 cmH₂O. This ensures that the cuff is in close contact with the tracheal wall. Release the pressure immediately by pressing the red release valve until the pointer reaches the green area. We recommend the intra cuff pressure to be above 22 cmH₂O due to the risk of aspiration or pneumonia and below 32 cmH₂O due to the risk of ischemia of the tracheal mucous membrane (see publication of J. Stauffer in Respiratory Care, July 1999). Check and monitor the cuff pressure continuously. Adjust pressure increase via the red release valve and pressure decrease via the inflation bulb.

EXTUBATION

Before extubation attach the inflation line to the vacuum valve and squeeze the inflation bulb until it does not return to its original position. In order to check for complete deflation, attach the inflation line to the Luer nozzle of the cuff pressure gauge. The needle of the Cuff Pressure Gauge has to show a negative pressure by moving counter-clockwise. Only then the cuff is completely deflated and the patient can be extubated.

WIPE DISINFECTION

Wipe disinfection using conventional surface disinfectants based on alcohol (e.g., Incides® N from Ecolab Deutschland GmbH) is suitable for disinfection. When choosing a product for disinfection, a disinfectant with appropriate ranges of action must be used: bactericidal, yeas-ticidal, tuberculocidal and virucidal. After the wipe disinfection, the devices must be inspected for visible contamination. If necessary, repeat the wipe disinfection. After the wipe disinfection, the devices must be checked in accordance with the section "Functional check".

CAUTION

- Do not immerse the product in liquids.
- The product cannot be automatic or manual cleaned respectively sterilised.
- Do not clean the glass (plastic material) with aggressive agents.

METROLOGICAL INSPECTION

The accuracy ± 2 cmH₂O of the Cuff Pressure Gauges must be checked every 24 months. If required measuring equipment is not available, the metrological control can be carried out by the manufacturer. In this case the Cuff Pressure Gauge has to be sent to the manufacturer.

If there is an indication that the Cuff Pressure Gauge cannot keep the accuracy of measurement within the 24 months, or the metrological characteristics of the Cuff Pressure Gauge have been affected, a metrological control has to be carried out immediately. This has to be carried out only by a person who meets the requirements for metrological controls. The results of the control, as well as the measured values and measurement procedures has to be documented.

After successful control, the Cuff Pressure Gauge has to be marked with a sign which indicates (clearly and traceable) the year of the next inspection and the person who carried out the metrological control.

In order to check the values, the Luer nozzle of the Cuff Pressure Gauge has to be connected to a calibrated pressure gauge (measurement inaccuracy max. ± 0.6 cmH₂O).

Afterwards inflate via the inflation bulb to 30, 60 and 90 cmH₂O and record the values. The values must be within the indicated tolerance. If this is not the case, the device must be returned to the service of the manufacturer. If the values are within the tolerance, the device can still be used.

Note for conversion of cmH₂O to SI-Unit Pascal:

$$1 \text{ cmH}_2\text{O} = 0,98067 \text{ hPa}$$

SERVICE

Prior to return of medical products for complaint/repair, they must be cleaned and disinfected to exclude any risk for the manufacturer's staff. For safety reasons, the manufacturer preserves the right to refuse soiled or contaminated products.

STORAGE AND TRANSPORT CONDITIONS

The shelf life of the Cuff Pressure Gauge is 8 years.

The shelf life of the connecting tube is 5 years.

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

MAINTENANCE

The medical product may only be serviced by persons who have the technical knowledge, qualification and the required means to properly carry out the service.

After service of the medical product, the essential constructive and functional characteristics

for safety and functionality have to be checked, insofar as they can be influenced by the service procedure.

To ease repair of the medical product, return it together with a detailed description of the defect.

DISPOSAL

Used or damaged products must be disposed in accordance with the applicable national and international regulations.

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